



# SIERRA SCOOTER

## MEDIUM

### Classic Meets Comfort for Everyday Adventures

This scooter combines classic design with smooth operation to ensure a comfortable ride. Sierra gives you confidence in your daily travels thanks to its powerful motor and extended travel range.





### Stable and Secure on Every Trip

4-wheel independent suspension provides a steady ride over varied terrain.



### Extended Daily Travel

High-capacity lead-acid batteries deliver up to 40km of range per charge for peace of mind when out and about.



### Designed for Easy Manoeuvring

A tight 1565mm turning radius makes navigating narrow paths, shop aisles, and small community spaces simple.



### Plush, Supportive Seating

Features a deeply padded, adjustable captain's seat and flexible tiller, seat slides, swivels and reclines for perfect positioning.



### Simple and Intuitive Controls

Straightforward analogue interface is easy for new and experienced users to operate.



### Smart Safety

Emergency stop button and automatic speed reduction while turning keeps you safe in unexpected situations.

Overall L x W x H	1290 x 625 x 1150mm
Seat Size D x W	450 x 480 - 530mm
Seat to Deck Height	440 - 490mm
Seat to Ground Height	585 - 635mm
Armrest Height	215mm
Min Turning Radius	1565mm
Ground Clearance	87mm
Max Climbing Angle	12°
Max Driving Speed	10km/h
Driving Range	40km
Tyres	290mm   Pneumatic
Battery	2 x 12V 50Ah (Lead Acid)
Product Weight (excl. battery)	67kg
Battery Weight	30kg
Power (Nominal   Peak)	700W   1400W
Motor	4 Pole
Controller	RHINO 2 120A
Brake System	Electro-Mechanical
Free Wheel Mode	Yes
Charger Size	24V 5A



SCT671261  
**Arctic Mist**



SCT671271  
**Desert Sand**



\*Scooter Frame 7 Years | Motors/Axles/Brakes 5 Years | Electronics 1 Year | Battery 6 Months

**7 YEARS LIMITED WARRANTY**

**SWL 136kg**

170226