



# AeroLite Carbon Fibre Powerchair

Experience the freedom and comfort of the AeroLite powerchair. Its light carbon fibre frame is easy to fold, lift, and carry making loading into a car or storing away effortless.

Whether you're adventuring locally or overseas, AeroLite's lithium-ion batteries are approved for air travel, allowing you to explore the world freely.



## SIDE-FOLDABLE DESIGN

Stands when folded making lifting easy



## TWO SIZES AVAILABLE



Seat Width 400mm  
Blue PWS654145



Seat Width 450mm  
Red PWS654150



### Intuitive Joystick

Simple navigation for stress-free steering and greater control during daily activities



### Tool-free Assembly and Disassembly

Quick and hassle-free set up so you can start your journey quicker



### Two Long-range Batteries

Travel with peace of mind with a 30 km battery range (15 km per battery)

IATA-approved for air travel\*\*




### Durable Puncture-proof Tyres

Solid tyres allows users to travel with confidence over different surfaces



### Flip-up Footplates

Makes transferring in and out of the powerchair quick and easy



### Optional Accessories Available

Cup Holder SCA669200  
Phone Holder SCA669210

Overall L x W x H 960 x 630 x 895mm (Red) | 960 x 580 x 895mm (Blue)

Folded L x W x H 960 x 315 x 800mm

Seat Depth 400mm

Seat Width 400mm (Blue) | 450mm (Red)

Anterior Seat Height 480mm

Posterior Seat Height 440mm

Backrest Height 420mm

Backrest Width 400mm (Blue) | 450mm (Red)

Armrest Height 235mm

Armrest Length 350mm

Max Obstacle Height 25mm

Min Turning Radius 900mm

Max Climbing Angle 6°

Max Driving Speed 6km/h

Driving Range 30km

Front/Rear Wheels 180mm/300mm

Battery (Lithium-ion) 1 x 24V 12Ah (plus 1 x 24V 12Ah spare)

Battery Charger 2.6A

Motor/Drive System 150W

Min Product Weight (No Battery or Foot Plate) 13.2kg

Total Weight 16.6kg

Battery Weight 1.45kg (each)

Footplate Weight 0.5kg



**2** YEAR WARRANTY  
Frame

**1** YEAR WARRANTY  
Components

\*Range of motion dimensions are accurate to a +/-10mm tolerance

\*\*The AeroLite's Lithium-ion battery is compliant with international regulations. Please check with airlines for specific compliance regarding UN3171.