


# MONET DUAL ACTION

## Dual Action

The Monet Dual Action offers independent headrest function along with recline and tilt. The heating function provides warmth and comfort during the colder months whilst Zero Gravity tilt-in-space positioning creates a feeling of weightlessness and promotes pressure redistribution.



Canyon Ocean  
Microsuede



## Features & Benefits

- **Independent Adjustable Headrest** improves comfort whilst watching TV or reading a book.
- **Integrated Lift and Recline** functions for easy sit-to-stand.
- **Extending Legrest** supports more of the lower legs for enhanced comfort.
- **Heating Function** within the backrest provides warmth and comfort during the colder months.
- **Zero Gravity** tilt-in-space positioning creates a feeling of weightlessness and promotes pressure redistribution.
- **Backlit Hand Control** with integrated USB charging point.
- **Lithium Ion Battery Backup** ensures chair can be used in the event of a power outage or if moved between rooms during the day.

**5** YEAR WARRANTY  
FRAME

**2** YEAR WARRANTY  
COMPONENTS\*

**S.W.L**  
**159KG**

MONET DUAL ACTION CHP227540	SEAT HEIGHT (mm)	SEAT DEPTH (mm)	SEAT WIDTH (mm)	BACKREST HEIGHT (mm)	OVERALL HEIGHT (mm)	OVERALL WIDTH (mm)	OVERALL DEPTH (mm)	DISTANCE FROM WALL (mm)
	495	535	560	570	1065	890	1065	405

\* Components: Handset, upholstery and electronics.