


# REMBRANDT QUATTRO

## Quattro

The Rembrandt Quattro offers a unique seating experience. With independent control of back, leg, lumbar and headrest functionality, the user has infinite positioning options with the press of a button. Genuine shiatsu massage coupled with optional heating functions gives users a truly unique and luxurious comfort experience. Zero Gravity tilt-in-space positioning creates a feeling of weightlessness and promotes pressure redistribution.



Canyon Walnut  
Microsuede











## Features & Benefits

- **Shiatsu Roller Massage** in backrest provides therapeutic and stress relief for the user.
- **Memory Function** allows for a preset favourite position to be set.
- **Independent Adjustable Headrest** improves comfort whilst watching TV or reading a book.
- **Integrated Lift and Recline** functions for easy sit-to-stand.
- **Extending Legrest** supports more of the lower legs for enhanced comfort.
- **Heating Function** within the seat, lumbar and backrest provides warmth and comfort during the colder months.
- **Cup Holder** built into the armrest for easy drink storage.
- **Zero Gravity** tilt-in-space positioning creates a feeling of weightlessness and promotes pressure redistribution.
- **Wireless Device Charging** built into the armrest for convenience.
- **Lithium Ion Battery Backup** ensures chair can be used in the event of a power outage or if moved between rooms during the day.
- **4 x Swivel Castors** allow the chair to be moved around with ease.
- **Backlit Hand Control** with integrated USB charging point.

**5** YEAR WARRANTY  
FRAME

**2** YEAR WARRANTY  
COMPONENTS\*

**S.W.L**  
**159KG**

REMBRANDT QUATTRO CHP227550	SEAT HEIGHT (mm)	SEAT DEPTH (mm)	SEAT WIDTH (mm)	BACKREST HEIGHT (mm)	OVERALL HEIGHT (mm)	OVERALL WIDTH (mm)	OVERALL DEPTH (mm)	DISTANCE FROM WALL (mm)
								
	495	500	560	600	1095	1015	1065	405

\* Components: Handset, upholstery and electronics.