


STRIDE WALKER

EVERYDAY, LIGHTWEIGHT MOBILITY

Side folding walking frame, ideal for storage and transport



Removable utility bag



Reduced overall footprint for easy walker access



Ergonomic handles

Light touch lock brakes




Push button handle height adjustment

Soft padded seat for improved sitting comfort

Aluminium frame, lightweight and easy to manoeuvre

Multi-terrain castors



ISO 11199-2:2005

12 MONTH WARRANTY

RED

PRODUCT CODE WAF750150

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	FOLDED WIDTH	CASTOR HEIGHT	OVERALL WEIGHT	WEIGHT CAPACITY
790-980mm	830mm	600mm	230mm	480mm	520mm	230mm	Front - 11 inch Back - 9 inch	9kg	135kg