

# Jazzy<sup>®</sup>

## SELECT 6

The Jazzy Select<sup>®</sup> 6 incorporates sporty design and performance features for a truly unique ride experience. Pride's patented Active-Trac<sup>®</sup> Suspension combined with six wheels on the ground and an innovative articulating beam allow for superior performance indoors and out, with added stability over uneven terrain and ramp transitions. Mid-wheel drive provides the tightest turning radiuses in their class and caster wheels absorb the shock of obstacles, offering outstanding handling and control.

*Jazzy Select<sup>®</sup> 6  
shown in Jazzy Red*



**Pride**<sup>®</sup>  
Mobility Products  
Australia Pty. Ltd.

(Tel) 03 8770 9600 • (Fax) 03 9703 2720 • [www.pridemobility.com.au](http://www.pridemobility.com.au)



## SPECIFICATIONS

<b>Weight Capacity</b>	136 kg (300 lbs.)
<b>Maximum Speed</b>	Up to 6.8 km/h (4.2 mph)
<b>Ground Clearance</b>	4.45 cm (1.75") at motors and at center frame
<b>Turning Radius</b>	54.6 cm (21.5")
<b>Overall Length</b>	90.2 cm (35.5")
<b>Overall Width</b>	59.7 cm (23.5")
<b>Drive Wheels</b>	25.4 cm (10")
<b>Front Wheels</b>	15.24 cm (6") solid casters
<b>Rear Wheels</b>	15.24 cm (6") solid casters
<b>Suspension</b>	Active-Trac® (ATX)
<b>Drivetrain</b>	Two-motor, in-line, mid-wheel drive
<b>Braking System</b>	Regenerative and electro-mechanical
<b>Standard Electronics</b>	40 amp, PG GC 3 controller
<b>Battery Charger</b>	Off-board, 3 amp
<b>Per-Charge Range</b>	Up to 20.8 km (12.92 miles)
<b>Battery Requirements</b>	(2) 12 volt, deep cycle, U-1
<b>Battery Weight</b>	11 kg (24.5 lbs.) each
<b>Base Weight</b>	45.4 kg (100 lbs.)
<b>Standard Seat Weight</b>	19.5 kg (43 lbs.) (comfort, high-back)
<b>Standard Seat</b>	46 x 46-20 cm (18" x 18-20")

## FEATURES

- 6" castors on front & rear for maximum stability
- In-line motor technology for increased efficiency, torque and range
- Articulating rear castor beam
- Easy access charger port and circuit breaker
- Height and depth-adjustable foot platform
- Built with ease of service in mind
- 15.24 cm (6") castors on front and rear for maximum stability

## OPTIONS & ACCESSORIES

- Cup holder
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane/crutch holder
- Saddle bag

## COLOUR CHOICE



Jazzy  
Red

## OPTIONAL SEATING

- Synergy® Seat
- Solid seat pan

**Pride**  
Mobility Products  
Australia Pty. Ltd.

20-24 Apollo Drive Hallam, Victoria 3803  
ACN 088 609 661 ABN 66 088 609 661

(Tel) 03 8770 9600 • (Fax) 03 9703 2720 • [www.pridemobility.com.au](http://www.pridemobility.com.au)



The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.

